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Finger Stretching Trainer

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Instruction Manual

(SKD)

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Contents

Product Introduction	02
User Guide	04
Safety Precautions	13
Safety Precaution Include	14
Troubleshooting	16
Technical Parameters	17
CustomerService	18

After service, Thank you for choosing the finger trainer.

Please read this manual carefully before using this product.

Please keep this manual properly for guidance when using the stretching trainer.

The finger training device has two modes: the passive mode and mirror training mode. Our device aims to maintain and expand the range of motion of finger joints, prevent joint contracture and muscle atrophy, greatly promote the reorganization and compensation of the nervous system, and help patients recover as soon as possible.

The recovery of hand dysfunction after stroke is a long-term and arduous process. Patients' conditions, rehabilitation courses, and curative effects are very different. Under normal circumstances, the recovery of hand function is slower than other parts, and long-term rehabilitation training is needed.

If you have any questions in use, please contact: info@ziglite.com

The manual version may upgrade due to changes in software version or technical specifications without prior notice.


Product Introduction

Scope Of Application

It is suitable for rehabilitation training in (families) institutions for patients with finger joint recovery after injury or limb dysfunction sequelae caused by brain injury, stroke, cerebral palsy, spinal cord injury, orthopedic surgery, multiple sclerosis, etc.


Product Composition

Trainer Host


Product Introduction


Training Gloves


Mirror Glove

Power Adapter

Pressure Handle

User Guide

Buttons, Screen Icons In Detail

1

default to training mode

Training/Mirror Mode Display

The left image shows the working status indication of the training mode;

The figure on the right shows the working status indication of mirror mode.

(All green lights indicate the current mode status)

2

default to 9th gear

Strength Level Display & Settings

Press this button to set the training strength level, and the corresponding function display above will flash on the digital tube. At this time, you can use the "+" and "-" keys to set the range from 1 to 9 gears.

User Guide

3

default to 1st gear

Finger Opening Time Level Display & Settings

Press this button to set the finger opening time level, and the corresponding function display above will flash on the digital tube. At this time, you can use the "+" and "-" keys to set the range from 1 to 9 gears.

4


default to 1st gear

Finger Contraction Time Level Display & Settings

Press this button to set the finger contraction time level, and the corresponding function display above will flash on the digital tube. At this time, you can use the "+" and "-" keys to set the range from 1 to 9 gears.

5

Training Time Display & Time Setting

This knob can be used to set the training time, with a range of 15-45 minutes and a graduation value of 5 minutes.

User Guide

Buttons, Screen Icons In Detail

6

Return Button

Press this key to return the machine to the default mode of initial startup in any mode.

7

Mirror Mode Button

Press this button to enter/exit the mirror training mode.

8

+ "Setting Button

Press this button to increase the settings for training strength, finger opening time, and finger contraction time.

9

- "Setting Button

Press this button to reduce the settings for training strength, finger opening time, finger contraction time, and training time.

10


Start / Pause button

Press this button to start and pause training settings.

User Guide

Buttons, Screen Icons In Detail

11

Training Finger Display&Selection

Press these 5 buttons to select training fingers. When the blue light is on, it indicates that the finger training has been selected. When the selected finger LED light is not on, the finger does not participate in the training plan.

12

Action Mode Switching

Press this button to switch the action mode:

- ① When the blue LED lights on the fingers are all on, it indicates that they are in the training mode of 5 fingers together;
- ② When the blue LED light on the finger lights up like a walking lantern, it indicates that it is in finger by finger training mode.

Note: The training mode controls the timing of fingers through program control, which first controls the four fingers outside the thumb to contract first, and the thumb will delay contraction by 2 seconds. The advantage is to avoid causing the thumb and index finger to press against each other when the five fingers contract at the same time, which affects the training effect.

User Guide

Use Guide

1 Place The Trainer Device Host

Place the training host on the desktop, not on the wet, soft cloth or unclean desktop.

2 Connect Air Pipe And DC Plug

When connecting the training gloves, it is necessary to clip the air pipe connector into the corresponding interface of the host.

3 Wear Training (Mirror) Gloves

The guardian shall help the patient wear gloves and pull the Velcro glove to the wrist to bond and tighten.

User Guide

Use Guide

4

Power On

Turn down the "power switch" button of the trainer host to the "ON" position to turn on the power, and the display screen will light up to enter.

5

Start Hand Training

Press the "cycle mode / mirror mode" button to start and pause training at any time.

6


Shut Off

After the training, turn the "power switch" button of the trainer host to the "OFF" position to power off.

User Guide

Glove Wiring Diagram


Glove Joint Insertion Effect

User Guide

Air Valve Switch Control

1 Digital tube display

2 Function display

3 Control buttons

1 Pressure display:
The set range of contraction pressure range is 0-100KPa \pm 10%;
Set the range of opening pressure range to 0-60KPa \pm 10%;
Graduation value: 1KPa.

2

 The finger contraction pressure is set and displayed, and the corresponding function indicator light above is green. At this time, the setting range can be 0-100KPa \pm 10% through the "+" and "-" keys.

 The finger opening pressure is displayed for setting, and the corresponding function indicator light above lights up green. At this time, the setting range is 0-60KPa \pm 10%, which can be set using the "+" and "-" keys.

User Guide

<p>3</p> <p>Pressure Setting Function Selection</p> <p> When the host is powered on and turned on, the default pressure value of the handle is the maximum value. Press this button to set the opening/retraction pressure.</p>	
<p>“+” Set Button</p> <p> Press this button to increase the training pressure setting.</p>	<p>“-” Set Button</p> <p> Press this button to reduce the training pressure setting.</p>

- Note:** 1. When the host is powered on and turned on, the handle will also be displayed synchronously. At this time, the default pressure value of the handle when turned on is the maximum value state;
2. When pressing "  " to enter the pressure setting mode, regardless of whether the setting has been modified or not, it will run with the set value displayed on the pressure handle! To exit and cancel the pressure handle operation, press the "  " button. To re-enter the pressure handle operation mode, press the "  " button of the pressure handle again;
3. When entering the operating pressure handle mode, the power level key on the host will be closed and cannot be operated.

Safety Precautions

[warning] Failure to comply with the following safety instructions may lead to severe accidents such as fire, electric shock, and casualties.

- Do not disassemble, repair, or modify the components of the equipment by yourself to avoid abnormalities in the treatment device.
- Do not store components in or near water heaters, microwave ovens, cooking equipment, or high-pressure containers to avoid fire caused by overheating.
- It cannot be sterilized or autoclaved.
- Do not continue to use it when the alarm prompt appears.
- It is suitable for hand circulation and mirror movement rehabilitation training. Please do not use it for other purposes to avoid accidents.
- Do not store the parts in a place that is too hot or cold to avoid accidents.
- The validity period of the product is five years, and the disposal of the product at the end of its service life should comply with the requirements of local laws and regulations.
- Please keep the power connector dry and water-free.
- Please do not plug in or unplug the power supply with wet hands to avoid electric shock or injury.

Working Environment

- 1) Ambient temperature: 10°C~35°C
- 2) Relative humidity: 30%~75%, excluding condensation
- 3) Atmospheric pressure range: 700hpa~1060hpa
- 4) Power supply: AC100~240V, 50/60Hz (trainer host power supply DC12V)

Storage And Transportation

- During transportation, please fix the trainer in a relatively stable place to prevent collision and damage to the trainer.
- Please prevent liquids or foreign objects from entering the trainer to avoid damaging internal electronic components and causing malfunctions.
- Please do not place it in places with strong sunlight or ultraviolet radiation for a long time to avoid causing discoloration of the shell.
- Do not place other items on the product when saving it.
- Do not fold the gloves forcefully to avoid damaging the soft-actuated actuator and causing malfunctions.
- Do not store the trainer in a humid environment or magnetic solid field interference unsuitable for working with electronic devices.

Safety Precaution Include

- In case of skin allergy and other discomfort while using the product, please stop using it and consult a doctor.
- In case of a sprain, trauma, bone injury of unknown condition, severe osteoporosis, untreated fracture, or people who need bracing after treatment, do not use this product without professional confirmation.
- For patients with diseases or wounds or being treated for them or with a previous medical history, please consult a doctor after becoming familiar with the following contents.

Contraindication

- Patients with severe blood flow disorder such as thoracic venous tumors and hemocolumn disease.
- Patients with ulcer, wound, purulent fruit, bleeding, scald, inflammation, etc.

Other Precautions

- 1) In case of abnormal experience or body feeling during use, stop using it and use it under the guidance of a doctor after seeing a doctor;
- 2) In case of redness, pruritus and rash during use, stop use and seek medical attention;
- 3) Do not use it together with other equipment or in the state of applying drugs;
- 4) Do not use it alone for unconscious people or young children;
- 5) Do not use this product for any purpose other than the purpose of use.

Before Use

- 1) Before wearing training gloves, please try to run the therapeutic device host for 1 minute to confirm that all parts function usually;
- 2) Patients with mandatory systolic spasm need to massage before using the therapeutic device;
- 3) Note: before training, please confirm that the non muscle tension of the hand to be trained is too high (above grade II)

Safety Precaution Include

In Use

In case of the following conditions during operation, please stop using immediately and disconnect the power supply.

- 1) The trainer host emits abnormal sound;
- 2) The surface temperature of the trainer host is too high;
- 3) Water ingress into the trainer host;
- 4) Users experience pain, swelling, or other discomfort symptoms;
- 5) If any spasms occur during the training process, immediately stop using it.

Usage Time

- 1) Suggest selecting appropriate training methods based on one's actual situation;
- 2) The duration of each use should not be too long. If there is soreness in the hands the next day after training, it is necessary to reduce the number and duration of daily training;
- 3) It is recommended to use it 2-4 times a day, each time for about 15 minutes, with a minimum interval of 45 minutes between uses;
- 4) If it is necessary to replace or extend the use time during use, please consult professional medical personnel;
- 5) Long term continuous use can easily lead to high internal temperature of the machine, which may cause phenomena such as sleep or asynchronous movement. After the machine is fully cooled, it will automatically return to normal.

Use Of Power Adapter

- 1) The power adapter must be plugged and unplugged by a patient's supervisor;
- 2) After using the trainer, please disconnect the power adapter from the 220V power supply;
- 3) Do not use a power adapter that is not equipped with this machine.

Troubleshooting

Problem	The Reason	Measures
Flip the "power switch" down to turn on the machine	Poor power connection	Check the power cord and power adapter, the connection between the power adapter and the host
Training gloves are not working	Gloves that are not properly attached can damage the host or glove structure	Check the physical status of the training gloves and the connection between the trainer host and the training gloves

Note: If you still can't use the above method, please contact us.

Care & Maintenance

How to clean ?

- Do not wash training gloves and mirrored gloves directly with water. Dip a soft cloth in water or neutral detergent, wring it out, and wipe the gloves.
- When dust or dirt appears on the surface of the trainer host, users can use medical alcohol with a concentration of 75% to wipe it. When wiping, please dip a small amount of alcohol into a dry cloth to avoid allowing alcohol to drip or flow into the product.
- Clean with clean and dry soft cloth, and dry in a cool and ventilated place.

Note

- The use of cleaning and disinfecting agents other than alcohol may damage the product, shorten the service life or cause safety risks.
- It is forbidden to use high-temperature drying and exposure to the sun for quick drying or disinfection.
- Do not use gasoline, thinner or similar organic solvents to clean any part of this product. Do not use radiation, steam, ethylene oxide and other methods for disinfection, otherwise the product will be damaged.

Technical Parameters


Product Specifications

Product Name	Finger Stretching Trainer	Waterproof Level	IPX0
Product Model	SKD	Force Setting	1-9
Security Type	I Type	Time Setting	15-45 Minutes
Power Supply	AC100-240V 50/60Hz	Running Noise	≤60dB
Size	251*204*111mm	Training Mode	Cycle Training/ Mirror Training/Finger Selection Training

Glove Size

Size	S	M	L	XL	XXL
Hand Length (mm)	140-160	160-180	165-190	190-220	220-230
Hand Width (mm)	60-80	70-90	80-100	100-120	120-130
Middle Finger Circumference (mm)	45-50	50-60	60-70	70-90	90-100

- The corresponding model can be selected according to the actual hand measurement of the user. If two or more measured values are at the critical point, it is recommended to select one size larger.
- Women generally choose s/m, men generally choose l/xl, and the larger hand is XXL.

Customer Service

Email : info@ziglite.com

Address : Room B3, 18/F Bonsun Industrial Building,364-366 Sha Tsui Road,
Tsuen Wan, Hong Kong

Packing List		
Number	Name	Quantity
1	Trainer Host	1 pcs
2	Training Gloves	1 pcs
3	Mirror Gloves	1 pcs
4	Pressure Control Handle	1 pcs
5	Power Adapter Power Cord	1 pcs
6	Product Manual	1 pcs
7	USB Data Cable	2 pcs

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Office Address : Room B3, 18/F Bonsun Industrial Building, 364-366 Sha Tsui Road, Tsuen
Wan, Hong Kong

Email : info@ziglite.com

Fax: +852 3544 0462

Office Tel : +852 3107 7500

Whatsapp : +852 6827 2010

NO. SKD02-CH-EN

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Ziglite Smart Health Care Products Co., Ltd.